

C9orf163-chromosome 9 open reading frame 163 Gene

Reference: PTXBC117152

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C9orf163; chromosome 9 open reading frame 163		
For research use only.			

GENE DESCRIPTION	chromosome 9 open reading frame 163
GENE ID	158055
ORIGIN SPECIES	Human
ACCESSIONS	BC117152
VECTOR	pDONR223
NCBI SYMBOL	C9orf163
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggggccactcactgcacacgtgctggcaggcagggaagagccggccttctctgctgctcttcagaggcaggcgcgtgtgtgggggtcccggccgatggcatcgccggaggctgagctcccagcaggcgtgaggtcctcctgggtgggagccaccctt
gtccccagctggggccggcgtggtgaggagggggtgatatccgtgccacgccagcaggcgcaggcgtgagggagagcttcagtcggcgtgacgtggcaccaggcccaatctctcagtaacattgaccttcaggagtcagattcctgacctgtttgaacagagttaggga
gcacgtggtaggggccagcccgcccgctgccccatctgtcttccagtggtcagggcgtctgcaccctccgggaaggaggtgtcattgtttgccctttccaaaggaggatgcagaggtgagctgcctcagggcgtggggctgcctttgccgctcgagttcaag
aaccgaagagcagagccagggtcgtgacctctaggtcctccgaaccctgtgtgccaccgggtgggcactaa