

C22orf27-chromosome 22 open reading frame 27 Gene

Reference: PTXBC117160

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C22orf27; NCRNA00325; MORC2 antisense RNA 1		
For research use only.			

GENE DESCRIPTION	chromosome 22 open reading frame 27
GENE ID	150291
ORIGIN SPECIES	Human
ACCESSIONS	BC117160
VECTOR	pDONR223
NCBI SYMBOL	C22orf27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgctgagatgaagagagagcaggatggcagaggaatcaccagctcaagaacccaagacttggacagctctccagcagaagggcagggtctctcagtgactactcacaacaggcctctcctagagatgtagtctgggaagggtccaaaaccaccagacaggcacatc
cagtggggaagacttccattccagtgtggcctgcatccctcactctgcttccagttctgctgacctcagaggagcagggtggccattgctgactctgggttagggagatgccagaacacacccaggttgaaccacaggttagtggtacaaagtcacagtccttgccttcagaa
aaccagaggcccaaatcacgagtgaggagagagttgggaaaaactgaaccagactctccaatggccattaaagagaaaatcagcaacgcttgcctttgctttgacctgagaagttaa