



Reference: PTXBC117165

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPA7; CPAD; CYP2A; CYP11A7; P450-IIA4; cytochrome P450 2A7; cytochrome P450 IIA4; cytochrome P450, family 2, subfamily A, polypeptide 7; cytochrome P450, subfamily IIA (phenobarbital-inducible), polypeptide 7; cytochrome P450IIA4; cytochrome P450 family 2		

For research use only.

GENE DESCRIPTION	cytochrome P450, family 2, subfamily A, polypeptide 7
GENE ID	1549
ORIGIN SPECIES	Human
ACCESSIONS	BC117165
VECTOR	pDONR223
NCBI SYMBOL	CYP2A7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr