



Reference: PTXBC117171

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rho guanine nucleotide exchange factor 10-like protein; Rho guanine nucleotide exchange factor (GEF) 10-like; Rho guanine nucleotide exchange factor 10 like		

For research use only.

GENE DESCRIPTION	Rho guanine nucleotide exchange factor (GEF) 10-like
GENE ID	55160
ORIGIN SPECIES	Human
ACCESSIONS	BC117171
VECTOR	pDONR223
NCBI SYMBOL	ARHGEF10L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]