



Reference: PTXBC117175

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SMILE; transmembrane and TPR repeat-containing protein 3; transmembrane and tetratricopeptide repeat containing 3		

For research use only.

GENE DESCRIPTION	transmembrane and tetratricopeptide repeat containing 3
GENE ID	160418
ORIGIN SPECIES	Human
ACCESSIONS	BC117175
VECTOR	pDONR223
NCBI SYMBOL	TMTC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]