

FN1-fibronectin 1 Gene

Reference: PTXBC117176

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CIG; ED-B; FINC; FNZ; GFND; GFND2; LETS; MSF; cold-insoluble globulin; migration-stimulating factor; fibronectin 1		

For research use only.

GENE DESCRIPTION	fibronectin 1
GENE ID	2335
ORIGIN SPECIES	Human
ACCESSIONS	BC117176
VECTOR	pDONR223
NCBI SYMBOL	FN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]