



Reference: PTXBC117188

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BTB/POZ domain-containing protein KCTD3; NY-REN-45; NY-REN-45 antigen; potassium channel tetramerisation domain containing 3; renal carcinoma antigen NY-REN-45; potassium channel tetramerization domain containing 3		

For research use only.

GENE DESCRIPTION	potassium channel tetramerisation domain containing 3
GENE ID	51133
ORIGIN SPECIES	Human
ACCESSIONS	BC117188
VECTOR	pDONR223
NCBI SYMBOL	KCTD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]