

KRTAP19-4-keratin associated protein 19-4 Gene

Reference: PTXBC117197

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP19.4; keratin-associated protein 19-4; keratin associated protein 19-4		
For research use only.			

GENE DESCRIPTION	keratin associated protein 19-4
GENE ID	337971
ORIGIN SPECIES	Human
ACCESSIONS	BC117197
VECTOR	pDONR223
NCBI SYMBOL	KRTAP19-4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctactatggcagctattacagaggccctgggctatggctgtggaggtttgtggcctaggtatggctatggctgtggatgtggcagcttccgcagactgggttatggctgtggctttggaggcaacggatatggctgctgccgcccatcatgctatggagataggattctcaattct
actgaaatcctaccctgaggacacgataatcagaggtcataagaagatcattcaatttaaccaaatttga