



Reference: PTXBC117200

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PWD; WC1; WND; copper-transporting ATPase 2; ATPase, Cu(2+)- transporting, beta polypeptide; ATPase, Cu++ transporting, beta polypeptide; Wilson disease-associated protein; copper pump 2; ATPase copper transporting beta		

For research use only.

GENE DESCRIPTION	ATPase, Cu++ transporting, beta polypeptide
GENE ID	540
ORIGIN SPECIES	Human
ACCESSIONS	BC117200
VECTOR	pDONR223
NCBI SYMBOL	ATP7B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]