



Reference: PTXBC117222

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SID-1; SID1; SID1 transmembrane family member 1		
For research use only.			

GENE DESCRIPTION	SID1 transmembrane family, member 1
GENE ID	54847
ORIGIN SPECIES	Human
ACCESSIONS	BC117222
VECTOR	pDONR223
NCBI SYMBOL	SIDT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr