

ITLN2-intelectin 2 Gene

Reference: PTXBC117225

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HL-2; HL2; intelectin-2; endothelial lectin HL-2; intelectin 2		

For research use only.

GENE DESCRIPTION	intelectin 2
GENE ID	142683
ORIGIN SPECIES	Human
ACCESSIONS	BC117225
VECTOR	pDONR223
NCBI SYMBOL	ITLN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgtccatgctgaggacaatgaccagactctgtctctgtttattctctctgtggccaccagtgggtgcagtcagcagccctctctgtgagatgctctcaggggaattcgaaacctgtgccttctcttttctccctgcctagaagctgcaagaaatcaaggacgtgccatagtcag
gtgatggcctgtattttctccgaccaagaatgtgtgtgtaccagacctctgtgacatgacttctgggggtggcgtggacctgtgtgcccagctgcacgagaatgacatggaagtgcacgtgggtgacgtgtgtccagtcagcagggaacaaagcagactaccag
aggggatggcaactggccaactacaacctttggtctgagaggcggccagagcgtgactacaagaacctgtgctactacgacatccaggccaaggacctgggcatctggcatgtgcccaacagtcctccatgcagcattggagaacagcgcctgtgaggtaccgca
ccaacactggtctctccagagactgggacataatctgttggcatctaccagaataccagtgaaatacagatcagggaatgttggaaatgacaatggcccagccatacctgtgtgtctatgactttgtgtgatgctaagaagactgcactttattactcaccgtatgtcaacgggaattgttg
caggattcgttcagttccgggtgtttaataacgagagagcagccaacgccctttgtgtgtgataaaagtactggctgtaacactgagcatcactgcacgtgtggaggagggtcttctccacagggcaaacccgtcagttggggacttctccgctttgactgggatggaactc
acgttaagagcagctgcagtcgggagataacggagcggctgtactctgttctatagatga