



Reference: PTXBC117228

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NTPase KAP family P-loop domain-containing protein 1; NTPase, KAP family P-loop domain containing 1		

For research use only.

GENE DESCRIPTION	NTPase, KAP family P-loop domain containing 1
GENE ID	284353
ORIGIN SPECIES	Human
ACCESSIONS	BC117228
VECTOR	pDONR223
NCBI SYMBOL	NKPD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]