

HOXB5-homeobox B5 Gene

Reference: PTXBC117247

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HHO.C10; HOX2; HOX2A; HU-1; Hox2.1; homeobox protein Hox-B5; homeo box 2A; homeo box B5; homeobox protein HHO.C10; homeobox protein Hox-2A; homeobox protein Hu-1; homeobox B5		

For research use only.

GENE DESCRIPTION	homeobox B5
GENE ID	3215
ORIGIN SPECIES	Human
ACCESSIONS	BC117247
VECTOR	pDONR223
NCBI SYMBOL	HOXB5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctgtactttgaaactcttctcggggcgttatccaaatggcccgactatcagttgctaaattatggcagtgagctctctgagcggctcttacaggatcccgctgcatgcacacggctcttacggctacaattacaatgggatgacctcagcgtcaaccgtcctcggcctc
ctccagccactttggggcgtgggcgagagctcgcgcgcttccccgcggcccccaggagccccgcttcaggaagcggcttcgagctgtccctgtctcggccgagtcctcctgcctgcaccaaaggcgacagccagggccaagccctctgtctgtccccctccgaccaggc
gacctcagcagctccagcgcgaatttcacgaaatagacgagcgcgcgtcctcggagcctgaggaagcggcgaagccagctaaagcagccccagcctagctcgggcgagccagccatggccacctccacagccgcgcggcgggagactccgaaatattccctg
gatgaggaagcttcacatgagccatgatatgacccggcgacgggaaaggccggaccgcgtataccgctaccagaccctggagctggaaaaggagttccactcaaccgctacctgacccggcgacggcgatcagatgccacgcactctgcctgtccgagcgcaga
tcaagatctggttcagaaccggcgcatgaagtggagaaggaacaacaattgaaaatgatgagcctggctacagctggcagcgccttcagccctga