



Reference: PTXBC117255

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	histone deacetylase complex subunit SAP130; 130 kDa Sin3-associated polypeptide; Sin3-associated polypeptide p130; Sin3A associated protein 130kDa; Sin3A associated protein 130		

For research use only.

GENE DESCRIPTION	Sin3A-associated protein, 130kDa
GENE ID	79595
ORIGIN SPECIES	Human
ACCESSIONS	BC117255
VECTOR	pDONR223
NCBI SYMBOL	SAP130
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]