

C11orf30-chromosome 11 open reading frame 30 Gene

Reference: PTXBC117265

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf30; GL002; BRCA2-interacting transcriptional repressor EMSY; protein EMSY; EMSY, BRCA2 interacting transcriptional repressor		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 30
GENE ID	56946
ORIGIN SPECIES	Human
ACCESSIONS	BC117265
VECTOR	pDONR223
NCBI SYMBOL	C11orf30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr