



Reference: PTXBC117281

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein CASC1; LAS1; PPP1R54; cancer susceptibility candidate gene 1 protein; lung adenoma susceptibility 1-like protein; protein phosphatase 1 regulatory subunit 54; cancer susceptibility candidate 1		

For research use only.

GENE DESCRIPTION	cancer susceptibility candidate 1
GENE ID	55259
ORIGIN SPECIES	Human
ACCESSIONS	BC117281
VECTOR	pDONR223
NCBI SYMBOL	CASC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]