

SCN3B-sodium channel, voltage-gated, type III, beta Gene

Reference: PTXBC117282

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATFB16; BRGDA7; HSA243396; SCNB3; sodium channel subunit beta-3; sodium channel, voltage-gated, type III, beta subunit; voltage-gated sodium channel beta-3 subunit; sodium voltage-gated channel beta subunit 3		

For research use only.

GENE DESCRIPTION	sodium channel, voltage-gated, type III, beta
GENE ID	55800
ORIGIN SPECIES	Human
ACCESSIONS	BC117282
VECTOR	pDONR223
NCBI SYMBOL	SCN3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgcctcaatagattgttccctggcttctctgtcttacttgggtcagtgctgtcttccctgtgtgtggaagtgcctcggagacggaggcctgcagggaaccccatgaagctgcgctgcctctctcatgaagagagaggagtgaggccaccacgggtggaa
tgggttctacagggccgagggcggtaaagatttcttatttacgagatcgggaatgcccaccaggaggtggagagccctttcaggggcgctgcagtggaaatggcagcaaggacctgcaggacgtgtccatcaactgtgtcctaacgtcactctgaacgactctggcctctacacctgcaatgt
gtccgggagttgagttgagggcgcacggccctttgtgaagacgacggcgtgatccccctaaagagtcaccgaggaaggctggagaggacttcactctgtgtgtcagaaatcatgatgtacatccttctgtgttctcactctgtgtgtgtcctcgaatgatgatattgctacagaaagg
tctcaaaagccgaagaggcagccaagaaaacgctgtgactaccttgcacccatctgagaacaaggagaactctgcgtaccagtggaggaaatag