

IFNW1-interferon, omega 1 Gene

Reference: PTXBC117290

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	interferon omega-1; IFN-omega 1, interferon omega-1; interferon alpha-II-1		
For research use only.			

GENE DESCRIPTION	interferon, omega 1
GENE ID	3467
ORIGIN SPECIES	Human
ACCESSIONS	BC117290
VECTOR	pDONR223
NCBI SYMBOL	IFNW1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctctgttccctctactggcagccctagtgtgaccagctatagccctgttggatctctgggctgtgtatgcctcagaacctatggcctacttagcaggaacaccttgggtcttgcacaaatgaggagaatctcccttcttctgtctcaaggacagaagagactcaggtcccc
caggagatggtaaaaggagccagttgcagaaggcccatgtcatgtctgtctccatgagatgctgcagcagatcttcagccctctccacacagagcgctcctctgtgcctggacatgacctcttagaccaactccacactggacttcacagcaactgcaacacctggagacctgctt
gctgcaggtagtgggagaagagaatctgctgggcaattagcagccctgcactgacctgaggaggtacttccagggaatccgtgtctacctgaaagagaagaatacagcgactgtgcctgggaagttgtcagaatggaatcatgaaatcctgtttctatcaacaacatgcaagaa
agactgagaagtaagatagagacctgggctcatcttga