

DSCR8-Down syndrome critical region gene 8 Gene

Reference: PTXBC117304

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf65; CT25.1a; CT25.1b; MMA-1; MMA-1a; MMA-1b; MMA1; MTAG2; Down syndrome critical region gene 8		

For research use only.

GENE DESCRIPTION	Down syndrome critical region gene 8
GENE ID	84677
ORIGIN SPECIES	Human
ACCESSIONS	BC117304
VECTOR	pDONR223
NCBI SYMBOL	DSCR8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggagcctggacccaactttgttactgtgagaaagggtcttcattcattcaagatggcatttgttaagcacctactgctggagtcagtggttcaatcacggatcactgcagcctccacctccagttcaagaaattctcatgtctcagcctcctgagcagctaggattacagacaaacctt
ggaaatcaagaaagtctggaatgatgaagctgttcattgccaagaccgaaagtctggtgcccagatgatgagtcattcagttcatgccgtga