



Reference: PTXBC117329

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GC1; GW112; OLM4; Olfd; UNQ362; ba209J19.1; hGC-1; hOLfd; olfactomedin-4; G-CSF-stimulated clone 1 protein; antiapoptotic protein GW112; olfactomedin 4		

For research use only.

GENE DESCRIPTION	olfactomedin 4
GENE ID	10562
ORIGIN SPECIES	Human
ACCESSIONS	BC117329
VECTOR	pDONR223
NCBI SYMBOL	OLFM4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]