

PLAC9-placenta-specific 9 Gene

Reference: PTXBC117332

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	placenta-specific protein 9; placenta specific 9		

For research use only.

GENE DESCRIPTION	placenta-specific 9
GENE ID	219348
ORIGIN SPECIES	Human
ACCESSIONS	BC117332
VECTOR	pDONR223
NCBI SYMBOL	PLAC9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccgcctgctgctgacccgactggccctgtccgcgcgcgggctctttggccgctgccgaaccttcagccctccgcgaggagactcagctcagagcacagcgtgtgacagacacatggctgtgcaacgccgtctagatgtcatggaggatgtagagaagacc
gtggatcaccctgggacagaggtgaaaggcctgctgggctgctggaggagctggcctggaacctgccccgggacccttcagccccgctccgaccttctcgagatggcttctga