

C4orf26-chromosome 4 open reading frame 26 Gene

Reference: PTXBC117342

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C4orf26; AI2A4; amelogenesis imperfecta type IIA4; chromosome 4 open reading frame 26		

For research use only.

GENE DESCRIPTION	chromosome 4 open reading frame 26
GENE ID	152816
ORIGIN SPECIES	Human
ACCESSIONS	BC117342
VECTOR	pDONR223
NCBI SYMBOL	C4orf26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcagacactgcttctctactggttactggtatgctggttggtgtaactgtggcagaaggacaagaagggtatttacgctcctggagattcacaataatgcggacgtaccgactgccagatctttacactcaccctccacctgccccgaggagtcggtcacaagggcc
cagcccatcacaagacaccaggtgtcccttccattttttccacgaaggccagaatccatttttaggtttccaacagacctttcctcctcaagggtgaaccaccgttttccattccagccattttattggccacaccgttaccttactttaggtattttcccagaagaagactccagagagga
agctcatctgaggaaagctga