



Reference: PTXBC117359

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Kv7.5; potassium voltage-gated channel subfamily KQT member 5; KQT-like 5; potassium channel protein; potassium channel subunit alpha KvLQT5; potassium channel, voltage gated KQT-like subfamily Q, member 5; voltage-gated potassium channel subunit Kv7.5; p		

For research use only.

GENE DESCRIPTION	potassium voltage-gated channel, KQT-like subfamily, member 5
GENE ID	56479
ORIGIN SPECIES	Human
ACCESSIONS	BC117359
VECTOR	pDONR223
NCBI SYMBOL	KCNQ5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]