

LOC200261-hypothetical protein LOC200261 Gene

Reference: PTXBC117361

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC200261
GENE ID	200261
ORIGIN SPECIES	Human
ACCESSIONS	BC117361
VECTOR	pDONR223
NCBI SYMBOL	LOC200261
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagaagccacctgccccacagggggccgttccctgcatctggcactcaccaggacagtcataaactggcctgggtgcttggggagagtgaggaggtctggggggcatccactggcctcagcctcagagcttctcatcagtaaccaggaggcaaagggccagctcgaga
ttagccaggccactgtatctgtgccaagacctacagctgtccctgcaggcgctgcatctctgcagcccgctatcacaatgggtgtcccttttgggggatgtga