

LOC145757-hypothetical LOC145757 Gene

Reference: PTXBC117365

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC145757
GENE ID	145757
ORIGIN SPECIES	Human
ACCESSIONS	BC117365
VECTOR	pDONR223
NCBI SYMBOL	LOC145757
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctgccccagagtcgattttcttggagcaacgcaaagcacagcttctgggagaagacacaatgcctcctattgtctctgccacttccctccagcctccagcattttctgctgtctcctgggagccactgccctgtgcaagagctgtttgctgatgccagccgcttccaaac
tcccgcgagaagtcgtttctgtccacagtcattgttccagaccagccgagccaactggggagggagagccagggccagggagagcagtaagagggccagcttgcaacacagggccccggcgctggagtcggctcatgggttctggggggcagctggtgccagccctgctagac
atgtccagggctggggcagagggcagggtggctctcatccacaacgaggtctccaggaagggctccagagggccaggtcctggaggaagacagaggggtgattatccagtga