

ADCY10-adenylate cyclase 10 (soluble) Gene

Reference: PTXBC117366

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HCA2; HEL-S-7a; SAC; SACI; Sacy; hsAC; adenylate cyclase type 10; 3',5'-cyclic AMP synthetase; AH-related protein; ATP pyrophosphate-lyase; epididymis secretory sperm binding protein Li 7a; germ cell soluble adenylyl cyclase; testicular soluble adenylyl c		

For research use only.

GENE DESCRIPTION	adenylate cyclase 10 (soluble)
GENE ID	55811
ORIGIN SPECIES	Human
ACCESSIONS	BC117366
VECTOR	pDONR223
NCBI SYMBOL	ADCY10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr