

DEFB118-defensin, beta 118 Gene

Reference: PTXBC117378

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf63; DEFB-18; ESP13.6; beta-defensin 118; defensin, beta 18; epididymal secretory protein 13.6; epididymus specific clone 42; defensin beta 118		

For research use only.

GENE DESCRIPTION	defensin, beta 118
GENE ID	117285
ORIGIN SPECIES	Human
ACCESSIONS	BC117378
VECTOR	pDONR223
NCBI SYMBOL	DEFB118
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaactctgctgctggctctctatgctgtgctcctacccaagtgtatccagcctatagtggtgaaaaaaatgctggaacagatcagggcactgcaggaaacaatgcaaagatggagaagcagtgaaagatacatgcaaaaatcttcgagcttgcctccatccaatgaagac
cacaggcgagttctctcgacatctccacacccttgagtgactcaacaccaggaattattgatgatatitaaacagtaaggttcacgacagactactttgaagtaagcagcaagaagatatggttgaagagcttgaggcgggaagggaactgagacctctcttccaaatgttcaccatagct
catga