

COMMD6-COMM domain containing 6 Gene

Reference: PTXBC117391

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Acrq; COMM domain-containing protein 6; Acrq embryonic lethality minimal region ortholog; COMM domain containing 6		

For research use only.

GENE DESCRIPTION	COMM domain containing 6
GENE ID	170622
ORIGIN SPECIES	Human
ACCESSIONS	BC117391
VECTOR	pDONR223
NCBI SYMBOL	COMMD6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcgtccagcgagccgctggatgctaaagtcacacaccagctgtagattttcagtggaactgggtatggctgtgagctcagacactgcagatctctaaagtatcctacgttcagtgatgctaaaagtgagcagatcattcaggccaagtaagaccaagtccttg
aaatgacgattccagatttcagaattctacagacagttcaaggaaattgctgcagttattgaaacgggtgga