

CSTL1-cystatin-like 1 Gene

Reference: PTXBC117396

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CTES1; RCET11; dJ322G13.4; cystatin-like 1; cystatin like 1		
For research use only.			

GENE DESCRIPTION	cystatin-like 1
GENE ID	128817
ORIGIN SPECIES	Human
ACCESSIONS	BC117396
VECTOR	pDONR223
NCBI SYMBOL	CSTL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggatcggatgctggagaaacccctgctgctgattgccctggctctgagccaaagctgggtcacttccaaagggtggagggttcagcagaagctcatgagcaagaagaacatgaattcaacactcaactcttcattcaatctacaacaatgccagcaacgacacctactta
tatcgagtcagaggctaatcgaagtcagatgcagctgacgacgggagtgagtagatagtcactgtgaagattggccggacaaatgcaagagggaatgacacgagcaattctctgccccctgcaaagcaagaagctgagaaagagttaatttgcgagctcttgatataccatgcc
ctggataaactatttcagctctggaacaattcctgtctggaggccgagcatgtgggcagaaacctcagatga