

C21orf49-chromosome 21 open reading frame 49 Gene

Reference: PTXBC117397

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 21 open reading frame 49		

For research use only.

GENE DESCRIPTION	chromosome 21 open reading frame 49
GENE ID	54067
ORIGIN SPECIES	Human
ACCESSIONS	BC117397
VECTOR	pDONR223
NCBI SYMBOL	C21orf49
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatcatgtcaggcagcttatgatccaattgtcccatggccttaattccacaacctcatcatctacccttggcttccagaaatcatgactctcaatgttgaagaatggcctactgcatggaaggtgctgatccttttagaggcagccgtggaagagaaaagatcagaagagaaa
aggatattggttgcggaacatgtggacaagaagtccccaagaattgtga