



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TGIF2LX-TGFB-induced factor homeobox 2-like, X-linked Gene

Reference: PTXBC117398

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein TGIF2LX; TGIFLX; TGF-beta-induced transcription factor 2-like protein; TGFB-induced factor 2-like protein, X-linked; TGFB-induced factor 2-like, X-linked; TGIF-like on the X; TGFB induced factor homeobox 2 like, X-linked		

For research use only.

GENE DESCRIPTION	TGFB-induced factor homeobox 2-like, X-linked
GENE ID	90316
ORIGIN SPECIES	Human
ACCESSIONS	BC117398
VECTOR	pDONR223
NCBI SYMBOL	TGIF2LX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagccgctcgcggacggcccgctgtagaccaaaagccgggtgaaaagacagccggcggaagaccaaaagccgacacctaatactgctcgagaaataacgcagatagcagcgagatcttctcctaccagagcaacaagaagagcgcaagggaacttgcca
gcccagtcggttaagatctccgcgactgcatgtataagcatcggttaaggcctaccctcagaaagaagagaagcaaatctgtgcagagaagaccatttcttctgtgcagatttctaactggtttatcaatctcgcagacgcattctcccggatagcttcaacgcgtgagaacgacccc
atcattggcccaaaaaggcgcaagatgcccatgcacccactcgcagagcaccgaggcgctgtgctggccgaagtcaggccgacgttgctcagacaatgtacaaagcctgccctctggcccttgcgaaggcgagatgtcagagagaagcaaccagatcgagtcgcccc
tagccagaagctaccggaaatagacacggcgaagaaaaaggtcaaggtttatcatcatcccgtctctcagaactgtgtctcagaggagcacgccacttcagcagcttctctgctagtctgatcgacgagtaaaagggtcctccgactgctgagctagagaagaagcaagagc
ctaactaag

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1