



Reference: PTXBC117403

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BIMP1; CARMA3; caspase recruitment domain-containing protein 10; Bcl10 binding protein and activator of NFkB; CARD-containing MAGUK 3 protein; CARD-containing MAGUK protein 3; carma 3; caspase recruitment domain family member 10		

For research use only.

GENE DESCRIPTION	caspase recruitment domain family, member 10
GENE ID	29775
ORIGIN SPECIES	Human
ACCESSIONS	BC117403
VECTOR	pDONR223
NCBI SYMBOL	CARD10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr