

C14orf105-chromosome 14 open reading frame 105 Gene

Reference: PTXBC117405

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C14orf105; chromosome 14 open reading frame 105		
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 105
GENE ID	55195
ORIGIN SPECIES	Human
ACCESSIONS	BC117405
VECTOR	pDONR223
NCBI SYMBOL	C14orf105
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcctgagtcacttaagactaccttagggatgcaaaagtagcactttgcaaaacaaagagtagagactccctcgctggcctgtggactttgcattcaatcagaatttgaagaaaagactcatattcactggcaagactgcaggaccagaataaagccttgaaggcagct
gccacctttacaagaaaactggatggaagatattctacagcatccagagacatgtatttgacatcccttggaacacagagaaacagattatttaaaaggcatcccccacaaagacttcaaaagcttgaacccattgacttgcacgagtaattacttcaggaagactcctgagccagca
gaagccaggacaatgcacaaagcaaaaggtactagaaaagaaatgcaaaactcgaatgtatactctgagaacagacaataatgcatagatgcaagtctggaatgatccgtaaaagacaagagcccaatggagttaaaagaaagcttcatggagaggcaagaattaagcaa
agtcgaaggaccataaagccaagaaaaccttcaagcaccacaaaggaatgatgacatgaccttcaacctgttgcctgatgaaatcttgaacagaggtcccgaaattcaagaatacagaattttgaaacatcagcagtgataaactgctgtccctggaaaattggcaaatggaa
acatgcttcattgaacaagagccaggacagctctctgtggacagttccagctctgactcagatgagcaggggaagatgagaagaagccagcagcactggtgagaaccaggacagagagaatccacttttcgatgagtttttgcataagaataa