



Reference: PTXBC117420

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLC-6; chloride transport protein 6; chloride channel, voltage-sensitive 6; chloride voltage-gated channel 6		

For research use only.

GENE DESCRIPTION	chloride channel 6
GENE ID	1185
ORIGIN SPECIES	Human
ACCESSIONS	BC117420
VECTOR	pDONR223
NCBI SYMBOL	CLCN6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]