

TAS2R44-taste receptor, type 2, member 44 Gene

Reference: PTXBC117421

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TAS2R44; T2R31; T2R44; T2R53; taste receptor type 2 member 31; taste receptor type 2 member 44; taste receptor type 2 member 53; taste receptor, type 2, member 31; taste 2 receptor member 31		

For research use only.

GENE DESCRIPTION	taste receptor, type 2, member 44
GENE ID	259290
ORIGIN SPECIES	Human
ACCESSIONS	BC117421
VECTOR	pDONR223
NCBI SYMBOL	TAS2R44
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaactttatacccatcattttccagtggtgtagtggttctattgttattggaaatttgcataatggcttcagcattgtaaatccattgagcgggtcaagagacaaaagatctctttgctgaccagattctcactgctctggcggtctccagagtggttgcctgggtattattataaatt
ggatfcaactgtgttaatacagctttttatagtgtagaagtaagaactactgcttataatgtctgggcagtaaccggccatttcagcaactggctgctactagcctcagcatattttatgtcgaagattgccaatctcccaacttattttctcacttaagaggagagttaaagtgctattctg
gtgatgctgttggggcctttactatttttgcctgtcaacttttggtaaaacatgaagagattgtacggacaaaagaatgaaggaaactgacttgaagatcaaatgaggagtgacgtgaccttccagatgcgactgtaaccacgctaggaaacttagtcccttcactctgaccctgcta
tggttttgcgttaactgttctctgtgtaaacatctcaagaagatgcagctccatggtaaggatctcaagatccagcaccaggtccacataaaagcttgcacaaactgtgactttttctctgttatgtgccgtttactttctgtccataatgatatcagtttggagtttgggagctctgaaaacaa
acctgtctcatgtctgcgaagctattagattcagctatccttcaatccaccatcatcctgatttggggaacaagaagctaaagcagactttcttctcagtttgcggcaagtgggtactgggtgaaaggagagaagccttcatctccatag