

TAS2R43-taste receptor, type 2, member 43 Gene

Reference: PTXBC117423

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T2R43; T2R52; taste receptor type 2 member 43; taste receptor type 2 member 52; taste receptor, type 2, member 43; taste 2 receptor member 43		

For research use only.

GENE DESCRIPTION	taste receptor, type 2, member 43
GENE ID	259289
ORIGIN SPECIES	Human
ACCESSIONS	BC117423
VECTOR	pDONR223
NCBI SYMBOL	TAS2R43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgataacttttctacccatcttttccagtcgtgtagtggttacatttgtattggaatttgcataatggctcatagcactggtaattccattgagtcgttcaagagacaaaagatctcttggcgaacaaatctcactgctctggcggtctccagagttggttgcctgggtattattataaact
ggtaattcaactgtgtgaatccagcttttaatagtgtagaagtaagaactactgcttataatcttggcagtgatcaacatttcagcaactggctgctactacccctcagcatatttttgcctcaagattgccaaattctcactttattttctcacttaagaggagagtttaagagtgctattctgtg
tgatgttgggggccttgcatttttggcctgtcatcttttggataaacatgaatgagattgtgagacaaaagaattgaaggaacatgacttgaagatcaaatgaagagtgcaatgtacttttcaataatgactgtaacatggtagcaaaccttagtaccctcactctgaccctactatcttt
tatgctgttaactgttcttgtgtaaacatctcaagaagatgcagctccgtggaaggaatctcaagatccagcagcaaggtccacataaaagcttgcacaaactgtgatctcttctcttgttatgtgccattttctgtccataatgatacagtttgagagtttggagagctggaaaacaacc
tgcttcatgttctgcaaaagctattagattcagctatccttcaatcccccattcatctgatttggggaacaagaagctaaagcagactttcttctcagtttttggcaaatgaggtactgggtgaaaggagagaagacttcatctccatag