

C4orf6-chromosome 4 open reading frame 6 Gene

Reference: PTXBC117444

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C4orf6; aC1; expressed in neuroblastoma; chromosome 4 open reading frame 6		

For research use only.

GENE DESCRIPTION	chromosome 4 open reading frame 6
GENE ID	10141
ORIGIN SPECIES	Human
ACCESSIONS	BC117444
VECTOR	pDONR223
NCBI SYMBOL	C4orf6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacacccagaagcaaatcacaagactcacaatcaagaaccaatttttacaattttttctgtcagttgaatttgggaaggaggaaacacgcaaaaattttaccttcttcaattggacactatggacggaaatccaggagagctgaccttggaaactgcagacactcagataa
aacagagccagaatgcttgcagctggacctttgacccaaccccgatgtgactgtctccgggtcactcaactgcagtgcacatga