



Reference: PTXBC117453

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G-protein coupled receptor 35; KYNA receptor; kynurenic acid receptor; G protein-coupled receptor 35		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 35
GENE ID	2859
ORIGIN SPECIES	Human
ACCESSIONS	BC117453
VECTOR	pDONR223
NCBI SYMBOL	GPR35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1