

NPTN-neuroplastin Gene

Reference: PTXBC117462

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GP55; GP65; SDFR1; SDR1; np55; np65; neuroplastin; SDR-1; stromal cell derived factor receptor 1; stromal cell-derived receptor 1		

For research use only.

GENE DESCRIPTION	neuroplastin
GENE ID	27020
ORIGIN SPECIES	Human
ACCESSIONS	BC117462
VECTOR	pDONR223
NCBI SYMBOL	NPTN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgggttcgtcgtccagcgccttgcctctgttgcgtgtctgtgctccctccagggccaggcgcctcagaacgagccaaggttgcaccagtgaagaggtcattattcagagagccctgttccctgtcaccctcagtgtaacctcacctccagctctcac
acccttacatacagctactggacaaagaatgggtggaactgagtgccactcgtagaatgccagcaacatggagtacaggatcaataagccgagagctgaggattcagcgaataccactgcgtatcactttgtcagcgtctctaaagcaaacgccaccattgaagtgaagccgct
cctgacatcactggccataaacggagtgaagaagaatgaaggcaggatgccactatgtattgcaagtcagttggctacccccaccagactggatatggcgcaagaaggagaacgggatgccatggacattgtcaatcctctgtccgcttctcatcatcaacaagaaaattaca
ctgagttgaacattgtgaacctgcagatcacggaagaccctggcagtgatgaatgtaatgccaccaacgccattggctccgctctgtgttcactgtcctcaggggtcggagccacctggccccctctggccttcttgggaattctggtgaaattatcatccttggtgatcattgtgtga
tgagaagagggaagagccagatgaggttctgacgatgaaccagctggaccaatgaaaaccaactctaccaacaatcacaagataaaaacttgcgccagagaacacaaattaa