



Reference: PTXBC117472

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CHDC1; NP81; leucine-rich repeat and calponin homology domain-containing protein 1; leucine-rich repeats and calponin homology (CH) domain containing 1; neuronal protein 81; leucine rich repeats and calponin homology domain containing 1		

For research use only.

GENE DESCRIPTION	leucine-rich repeats and calponin homology (CH) domain containing 1
GENE ID	23143
ORIGIN SPECIES	Human
ACCESSIONS	BC117472
VECTOR	pDONR223
NCBI SYMBOL	LRCH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]