

IL29-interleukin 29 (interferon, lambda 1) Gene

Reference: PTXBC117482

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL29; IL-29; interferon lambda-1; IFN-lambda-1; cytokine Zcyto21; interleukin 29 (interferon, lambda 1); interleukin-29; interferon lambda 1		

For research use only.

GENE DESCRIPTION	interleukin 29 (interferon, lambda 1)
GENE ID	282618
ORIGIN SPECIES	Human
ACCESSIONS	BC117482
VECTOR	pDONR223
NCBI SYMBOL	IL29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcagctggaccgtggctgctgacttggctgtaggtggccgtggcagggccctgtcccaactccaagccaccacaactgggaaggctgccacattggcaggttcaatctctgtcaccacaggagctagcgagcttcaagaaggccaggagcgccttgaagagtca
ctcaagctgaaaaactggagttgcagctctctgtcttccccgggaattgggacctgaggcttccagggtgaggagcgcctgtggccttggaggctgagctggccctgacgtgaaggtcctggaggccgctgctggccagccctggaggacgtcctagaccagcccttcacac
cctgcaccacatctctccagctccaggccgtatccagctcagccacagcaggggccagggccggggccgctccaccactggctgcaccggctccaggagggcccaaaaaggagtcgctgctgctggaggcatctgtcacctcaaccttccgctctcacgcga
gacctcaaatatgtggccgatgggaacctgtgtctgagaacgtcaacccacctgagtcacctga