

NANOS2-nanos homolog 2 (Drosophila) Gene

Reference: PTXBC117484

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NOS2; ZC2HC12B; nanos homolog 2; NOS-2; nanos C2HC-type zinc finger 2		
For research use only.			

GENE DESCRIPTION	nanos homolog 2 (Drosophila)
GENE ID	339345
ORIGIN SPECIES	Human
ACCESSIONS	BC117484
VECTOR	pDONR223
NCBI SYMBOL	NANOS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagctgccaccctcgacatgtggaaggactactcaactgagccaggtggtgtggcgctgacgcaagtcggggctcaaaggctggagacccaagagattgaggagccaagtccgggctccgtggggcaggatcagggctggggcgccaggggccaacggggc
ctgggggacctgtgcaactctgcaagcacaacggggagtcggccagctctactctcacaccagctgaagacacgggatggcgtggtgtgtgtccatcctgaggcactacgtgtgtccgtgtgtcggggccaccgggtgaccaggccatagctcaagtactgccgcttaacgg
tgccagcagtcctctaccgccgagcggggcgcaactcgccgggacgcaggggtcaagcgctga