

NEUROG3-neurogenin 3 Gene

Reference: PTXBC117488

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Atoh5; Math4B; NGN-3; bHLHa7; ngn3; neurogenin-3; class A basic helix-loop-helix protein 7; protein atonal homolog 5		

For research use only.

GENE DESCRIPTION	neurogenin 3
GENE ID	50674
ORIGIN SPECIES	Human
ACCESSIONS	BC117488
VECTOR	pDONR223
NCBI SYMBOL	NEUROG3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgctcaacctcggtgcccactgtccaagtgaacctgagacggagcggctcctccagagcctcggaagacgaagtgaacctgccccacgtccgccccgcccagccactgcacacgggggaactgcgagaggcggaaggaggagctgccgagggggcccg
gaggaaagctccgggcacggcgccggggagcagccggcctaagagcgaagtggcactgagcaagcagcgacggagtcggcgaaagaaggccaacgacgcgcgaatcgaatgcacaacctcaactcggcactggagccctgcgcgggtgtctgtccaccttcccagac
gacgcgaagctcaccagatcagacgctgcgttcgcccacaactacatctggcgctgactcaaacgctgcgcatagcgaccacagctgtacgcgctggagccgcccggcgccgactgcgggagctgggcagccagggcggttccccggggactgggggtccctctact
ccccagctccccaggtggcagcctgagtcgccgcgtgctggaggagcgaccgggctgctgggggcccacctttccgctgttgagccagggcagctggtttctcagatttctgtga