

BATF3-basic leucine zipper transcription factor, ATF-like 3 Gene

Reference: PTXBC117489

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	JDP1; JUNDM1; SNFT; basic leucine zipper transcriptional factor ATF-like 3; 21 kDa small nuclear factor isolated from T-cells; 21-kD small nuclear factor isolated from T cells; B-ATF-3; Jun dimerization protein 1; Jun dimerization protein p21 SNFT; basic 1		

For research use only.

GENE DESCRIPTION	basic leucine zipper transcription factor, ATF-like 3
GENE ID	55509
ORIGIN SPECIES	Human
ACCESSIONS	BC117489
VECTOR	pDONR223
NCBI SYMBOL	BATF3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgcaaggctccggcgccggcagcgtcctgcagaggagcgtcgcggcgccgggaaccagccgagccgagccgagcagcagagccctgaggatgatgacaggaggtccgaaggagagaaaaaacccagttgctgctcagagaagtcggaagaagcagacc
cagaaggctgacaagctccatgaggaatatgagagcctggagcaagaaaacacatgctgcggagagagatcgggaagctgacagaggagctgaagcaccctgacagaggcactgaaggagcagcagagaagatgtgcccgctgctctgcctatgaactttgtccagtgccctc
ccggccggaccctgtggccggctgcttgcctccgatga