

LOC134466-hypothetical protein LOC134466 Gene

Reference: PTXBC117490

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC134466
GENE ID	134466
ORIGIN SPECIES	Human
ACCESSIONS	BC117490
VECTOR	pDONR223
NCBI SYMBOL	LOC134466
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcaggtctcagggtccagtgatgaggtgtggtgtggtttcaccaggaggagtgagcagcaactggactatgctcagaggaccctgtacagggatgtgatgctggagatctatagccactgggtctcaatgggatatccagttccaaccagatgtcatctcaagttggaac
aaggagaagagccatggatcataaagagacacataccaattggatctatccagacagagagatgagacttgacacccctcaactggatataattagagatgtttcttcataaaggagacactggaaagtattacaaggggtcattcattgtactcattttaaaagtctggcaaggtgatgac
cagctggagagagatcaggaaaactga