

RP5-1022P6.6-hypothetical LOC149837 Gene

Reference: PTXBC117492

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC149837
GENE ID	149837
ORIGIN SPECIES	Human
ACCESSIONS	BC117492
VECTOR	pDONR223
NCBI SYMBOL	RP5-1022P6.6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgcctttgcactgtccatgctctcagccaccacctcctccctcgcattgcagtggctaactggagacgaagaccaagaccacgcctttctccagtacctctcaaatcagcacagacaaggacaaggcctcaatcccaactgtgaagatggaccctggccacatgg
gatggtcagacacgcctgccagctatctgcaggcgaagaggtcagaagaggttttaggggcctgaaggacatctgtctccatgtccatagcaggtatttctgtccatga