

COL11A1-collagen, type XI, alpha 1 Gene

Reference: PTXBC117697

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CO11A1; COLL6; STL2; collagen alpha-1(XI) chain; collagen XI, alpha-1 polypeptide; collagen, type XI, alpha 1; collagen type XI alpha 1 chain		

For research use only.

GENE DESCRIPTION	collagen, type XI, alpha 1
GENE ID	1301
ORIGIN SPECIES	Human
ACCESSIONS	BC117697
VECTOR	pDONR223
NCBI SYMBOL	COL11A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr