

NWD1-NACHT and WD repeat domain containing 1 Gene

Reference: PTXBC117698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NACHT domain- and WD repeat-containing protein 1; NACHT and WD repeat domain-containing protein 1; NACHT and WD repeat domain containing 1		

For research use only.

GENE DESCRIPTION	NACHT and WD repeat domain containing 1
GENE ID	284434
ORIGIN SPECIES	Human
ACCESSIONS	BC117698
VECTOR	pDONR223
NCBI SYMBOL	NWD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1