



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C22orf34-chromosome 22 open reading frame 34 Gene

Reference: PTXBC118575

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

uncharacterized protein C22orf34; chromosome 22 open reading frame 34

For research use only.

GENE DESCRIPTION	chromosome 22 open reading frame 34
GENE ID	348645
ORIGIN SPECIES	Human
ACCESSIONS	BC118575
VECTOR	pDONR223
NCBI SYMBOL	C22orf34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]