

RP11-738I14.8-FLJ46082 protein Gene

Reference: PTXBC118593

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf171; cilia- and flagella-associated protein 77; cilia and flagella associated protein 77		

For research use only.

GENE DESCRIPTION	FLJ46082 protein
GENE ID	389799
ORIGIN SPECIES	Human
ACCESSIONS	BC118593
VECTOR	pDONR223
NCBI SYMBOL	RP11-738I14.8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagagggcaggagctccggcccgacacgcgatggaggaagcagcagcctgtgcggccacggcagcaggtctgcccccccgccggcccgctgaccgtggcgacatccgtccggcatggagaacgagcggctgggggtcgtcgggactcatgtt
cagaacccctctcatcgtaaggctgaactcggcaagccccgggaaagaagctacagctctgcccggcattatatttaattatggactctacatccgagggctgtgacggagagctccctgaagccatggagcgtggaacgtgtcaagcagcagccacctgccccacgagctgacccg
gaattatctgcaatgaaccgcccggcggtgaaagccggcctgtgactgcccgggagaactgtctacgtcagctcaatgacatcccatcagtgaccagagatgacggcgcatgaagaaagagccgccccctctccccaacatgacattgggatccggcagcgcttcca
caccctcttgatctgtgcagcaccggctacgtgcagctgtgggtacaggaacaaaggccaccagaaagccatcaaaactggagaagaagcagaaggtgtgtccttgggaagctgtatgagaccgggagcagtcagctgagggaagtacaagccgcccgtgaagctggacacccctg
gcacatgcctcacttcagaaggtggcgccaccttgatacgttccccacggagccgatccagagagcattaaaagccaccgggaagagtggtgccgtgcccaggggacccctgcggtggcgaactacacccacctag